NISSEI BLOOD PRESSURE MONITOR DSK-1011

INSTRUCTIONS

This manual is intended to assist you in the safe and efficient operation of BLOOD PRESSURE MONITOR DSK-1011. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

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

PL

CONTENTS

PART NAMES AND PRODUCT COMPONENTS	2
GENERAL INFORMATION	2
PRECAUTIONS AND CONTRAINDICATIONS	5
TIPS ON TAKING YOUR BLOOD PRESSURE	6
PREPARATIONS BEFORE USE	6
MAKING A MEASUREMENT	8
VIEWING THE SAVED READINGS (MEMORY FUNCTION)	10
TROUBLESHOOTING	11
TECHNICAL DESCRIPTION	13
WARRANTY	80

PART NAMES AND PRODUCT COMPONENTS

GENERAL INFORMATION

INDICATIONS FOR USE

This product is intended for noninvasive measurement of systolic and diastolic blood pressure, determination of pulse rate and calculation of pulse pressure in adults in a home healthcare environment. The product is not designed for neonatal use. Please consult with your doctor or physician to use this product to take blood pressure of child or person in pregnancy or under pre-eclamptic condition.

METHOD OF MEASUREMENT

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The cuff is connected to the main unit and wrapped around the arm. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and

displayed on the LCD as a digital value.

MEASUREMENT ON INFLATION

Initially the cuff is rapidly inflated to the prescribed value at the beginning of a measurement, during which the inflation mark is displayed. After the pressure reaches this point, the inflation mark disappears and the cuff is inflated more gradually, during which pulse wave is detected. As pulse wave is detected, heart mark flashes with synchronizing beep. Measurement ends when systolic blood pressure is determined and air is exhausted from the cuff.

PERSONAL MODE

The inflation is adjusted based on the latest three readings in the selected bank, either [1] or [2], to shorten the measurement time and to improve comfort.



For example, high diastolic blood pressure values from past three measurements will lead to inflation to higher pressure value than the default value. In the event that this value is too high for your current blood pressure, the cuff will rapidly deflate and measurement will restart.

To use the personal mode effectively, chose the memory bank before starting a measurement. This monitor remembers the last memory bank that was displayed. To choose the memory bank before starting a measurement, display a reading saved in the desired memory bank, turn off the monitor once and then start measurement.

WHO CLASSIFICATION SYMBOL

Measured blood pressure values are classified against WHO guideline (1999).

The higher class applies when systolic and diastolic pressures fall into different classes.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own based on the classification indication only.

DISPLAY	WHO CLASSIFICATION	SYS (mmHg)	DIA (mmHg)
	Hypertension (severe)	>180	>110
	Hypertension (moderate)	160-179	100-109
	Hypertension (mild)	140-159	90-99
	High Normal	130-139	85-89
	Normal	120-129	80-84
	Optimal	<120	<80

WHO: World Health Organization

PULSE PRESSURE DISPLAY

This product calculates and displays pulse pressure.

While systolic blood pressure keeps increasing with age, diastolic blood pressure tends to start decreasing at around age 50. Pulse pressure is value subtracting diastolic value from systolic value and therefore pulse pressure tends to increase with age. Pulse pressure has been considered to be related to stiffness of artery and has been studied as one of circulatory risk factors. It is sometimes said that 45 mmHg of pulse pressure is normal value. Pulse pressure only does not represent arterial sclerosis. However, observation of pulse pressure in long term will be important.


BODY MOTION SYMBOL

Blood pressure value taken while moving cannot be said to be the correct value because body movement can affect blood pressure.

This product analyzes pulse wave and displays [] when body motion is detected. [] indicates the results might be affected by body movement.

IRREGULAR PULSE RHYTHM SYMBOL

Pulse rhythm can be disturbed from talking, moving or arrhythmias. This product displays [] when the difference of shortest and longest time intervals between pulse waves is over 25%, indicating irregular pulse rhythm.

An example of irregular pulse rhythm compared to regular pulse rhythm

MEMORY FUNCTION

The measured values are automatically saved for later review in either one of two memory banks, [1] or [2]. These two memory banks can be used to save readings of two persons separately or to save morning readings and evening readings separately.

Each bank can save up to 60 readings. When the number of stored readings reaches 60, the oldest reading will be deleted to record new reading.

Saved readings are recalled with measurement date and time when clock function is activated. *The clock needs to be activated before measurement for date and time to be saved with the measured valued.

[Err] results are not saved.

CARE AND MAINTENANCE

Keep the product clean. Inspect its cleanliness after use. To clean, use only a soft dry cloth. Do not use gasoline, paint thinner, or other strong solvents. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, use a synthetic detergent and gently rub the surface. Air dry thoroughly. Make sure fluid never gets in AIR HOSE. Do not machine wash or scrub it.

When storing the product, do not place heavy objects on it. Do not fold or bend the cuff forcibly. To disconnect the cuff from the main unit, do not pull AIR HOSE but hold AIR PLUG and unplug it. Do not coil AIR HOSE too tightly. When the product has been stored

at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using. Remove the batteries if the product is to be stored for an extended period of time. Keep the batteries out of reach of children.

We suggest that you have your blood pressure monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

PRODUCT SPECIFICATIONS

Model	: DSK-1011
Operating Principle	: Oscillometric method
Indicator	: 15 digits liquid crystal display
Pressure Indicating Range	: 3 to 300 mmHg (cuff pressure)
Measuring Range	: 50 to 250 mmHg (systolic), 40 to 180 mmHg (diastolic), 40 to 160 bpm (pulse rate)
Accuracy*	: ± 3 mmHg (cuff pressure), $\pm 5\%$ of reading (pulse rate)
Inflation	: Automatic inflation
Exhaust	: Automatic exhaust valve
Power Supply	: Four 1.5 volt LR6 (AA alkaline) batteries or AC adaptor ADP-W5 series
Power Consumption	: 4W (max.)
Electric Rating	: DC6V/4W, with four LR6 batteries Adaptor: AC100-240V, 50-60Hz, 0.12A, Unit: DC6V, $=$ 500mA, with designated AC adaptor
Memory	: 2 banks, each stores 60 readings
Operating Condition	: +10°C to +40°C, 15% to 85% RH (noncondensing)
Transportation/Storage Condition	: -20°C to +60°C, 10% to 95% RH (noncondensing)
Cuff	: Coverage arm circumference: 22 to 32 cm, applied part; BF
Main Unit	: Weight: Approx. 250 g (without batteries) Size: Approx. 115 x 115 x 65.9 mm (W x D x H)
Protection class IP	: IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.

Protection against electric shock

: Internally powered equipment/Class II equipment, Type BF applied part

Mode of operation

: Continuous operation

Classification

: Internally powered equipment/Class II

Key to symbols

 : Type BF applied part

 : Class II equipment

 : Refer to instruction manual/booklet.

 : Keep dry

 : The used electrical and electronic products are not household

waste. Follow your national/local recycling rules to dispose of them properly. In the EU countries, please refer to waste management symbol(s) marked on the package or the instrument.

This device complies with EN1060-1:1995+A2:2009 Non-invasive sphygmomanometers Part 1: General requirements, EN1060-3:1997+A2:2009 Non-invasive sphygmomanometers Part 3: Supplementary requirements for electro-mechanical blood pressure measuring system and EMC (IEC60601-1-2:2007).

*Accuracy is guaranteed with the measured values that are within the measuring range. Specifications are subject to change without notice due to improvements in performance.

P R E C A U T I O N S A N D C O N T R A I N D I C A T I O N S

Do not use this product without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own regarding measurement results.

Do not use this product along with implantable and wearable medical electrical equipment such as pacemaker, defibrillator, or electrocardiographic monitor. This product is not also intended to be used with HF surgical equipment.

Do not use this product in an explosive environment such as near flammable anesthetics or inside oxygen chamber.

The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

Do not use cuffs or accessories other than those specified by the manufacturer. Otherwise, correct measurement readings cannot be obtained.

Use of this product in areas near mobile phones, microwave ovens or other

devices with strong electromagnetic field may cause malfunctions.

Do not apply the cuff over wounded arm, arm under an intravascular access or therapy or an arterio-venous shunt, or arm on the side of a mastectomy. Otherwise injury may be resulted.

Make sure that inflation of the cuff is not causing prolonged impairment of blood circulation. Also, be cautious about temporary loss of the functions of any other medical equipment if any monitoring equipment is used on the same limb with the blood pressure measuring cuff.

To avoid harmful injury due to interfered blood flow from cuff inflation,

- Make sure that AIR HOSE is not kinking before measurement. Otherwise, cuff inflation may not be conducted properly and prolonged, and
- Do not make measurements repeatedly.

To avoid any possibility of accidental strangulation, keep this product away from children and do not drape AIR HOSE around your neck.

Because the product includes precision parts, avoid extreme temperature variations, humidity, shock, dust, and direct sunlight. Do not drop or strike the

product. Make sure not to expose it to moisture. This product is not water resistant.

Do not press the display and operation keys or place the monitor with display face down.

Do not take out batteries or unplug the AC adaptor when the monitor is turned

on. Make sure to switch off the monitor before removing batteries or AC adaptor.

Do not touch the output plug of AC adaptor during measurement.

Do not disassemble or modify the product.

Do not inflate the cuff when it is not wrapped around your arm.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

Blood pressure varies all the time, influenced by mental and physical factors, and is never constant.

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement • Conversation • Mental Tension • Nervousness
- Emotions • Eating • Drinking Alcohol • Smoking
- Recent Urination or Bowel Movement • Room temperature
- Changes in the surroundings such as movement or noise, etc.

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.

This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

Let a qualified physician interpret your blood pressure readings.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Measure blood pressure using the correct posture and do not move or speak during measurement.

Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure at the same time every day.

The ambient temperature should be approximately 20°C when you take your blood pressure.

Reading may vary slightly, depending on posture, whether sitting or lying.

PREPARATION BEFORE USE

Install batteries and activate the clock before use. Measurement can be made with clock being inactivated. However, activate the clock to save measurement date and time with

readings for better review. The clock is inactivated as the batteries are removed.

The product can be also operated with AC adaptor. See page 8.

1 INSTALLING BATTERIES

Open the battery compartment cover, pressing the projection between [PUSH] and [█] at the bottom of the main unit.

Install four "AA" type batteries into the compartment.

Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment.

Batteries can be easily installed or removed as their (-) ends are pushed against the spring.

You may use nickel hydride rechargeable batteries.



NOTE: These batteries cannot be charged with this product. Use the designated battery charger. Also, for details regarding use of rechargeable batteries, follow the directions accompanying them.

You will hear a beep as the batteries are installed.

NOTE: The beep will not sound if AC adaptor is connected.

Close the battery compartment cover.

Do not force the battery cover into position. Insert two tabs to the main unit first, Fig.1). Then close the cover, pressing [█] section Fig.2). Close the cover completely Fig.3).

Battery symbol [] will flash when battery power is low. Measurement can not be made while the symbol is displayed. Replace the batteries. All four batteries need to be replaced with new ones; do not mix new and old batteries or different batteries.

The symbol might appear only during measurement although you can review readings in memory. This is because more battery power is required for conducting a measurement, e.g. inflating the cuff, than displaying the readings in memory. Please have the batteries replaced.

The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.


2 | ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

Touch "SET" KEY and keep your finger touched to the key until "20 11" starts flashing.

Clock is set in the order of year, month, day, hour and minute.


The flashing number increases with MEMORY KEY [1] and decreases with KEY [2]. The number will be fast-forwarded if you keep your finger on the key. Touching "SET" key will fix the number and the next article will flash. Touching "START/STOP" key will terminate the setting.

As the date and time are set, clock is activated. Display of the clock while the device is turned off indicates that the clock is activated.

OPERATION WITH AC ADAPTOR

AC CONNECTOR is located at the side of the monitor.

Use the designated AC adaptor only.

It is suggested to take out the batteries when using the adapter. However, the clock is inactivated as the AC adaptor is unplugged and therefore you may insert the batteries to

keep the clock activated. In such case, please make a few measurement without plugging the adaptor to avoid battery leak.

Do not connect the AC adaptor to any power socket where accessibility is not satisfying.

Protection: This device is double insulated and protected against short circuit and overload by a primary thermal fuse.


Enclosures and Protective Covers: Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test).

MAKING A MEASUREMENT


Insert AIR PLUG into AIR CONNECTOR before starting a blood pressure measurement.

The cuff is applicable to arm with circumference between 22 and 32 cm. Take your arm size before measurement.

Choose your memory bank before measurement to use the personal mode effectively.


Inaccurate readings may result if cuff is wrapped over heavy clothing or tight shirt that may restrict circulation in your upper arm or if sleeve is rolled up.


1 SIT ON A CHAIR AND APPLY CUFF.


Find a chair and a table so that you can be comfortably seated with your feet flat on the floor and your back and arm supported.

Put the cuff on your left arm with AIR HOSE positioned toward your hand.

Attach the cuff with hook and loop fastener with the lower edge of the cuff approximately 2 to 3 cm above the inner elbow.


The cuff should fit snugly, leaving enough space for two fingers to be put between your arm and the cuff.

If the cuff is wrapped more tightly or loosely than this, inaccurate blood pressure readings may result.

Adjust the position of the cuff so that AIR HOSE is over the inner part of your arm over the brachial artery.

Press the surface of the cuff to make sure that the hook & loop fastens securely.

2 Let your forearm rest naturally on the table and keep the cuff at the height of your heart during measurement.


If the cuff is lower (higher) than the heart, the measured reading tends to become larger (smaller).

Take deep breaths and relax.

Do not move or speak while taking your blood pressure.

Do not cross your legs during measurement.


Do not hold your breath.

3 TURN ON THE MONITOR.

Touch "START/STOP" KEY.

Air is exhausted from the cuff.


DEFLATION SYMBOL

Automatic inflation starts.

Touch "START/STOP" KEY again to cancel measurement.


The monitor will exhaust air from the cuff and turn off.

INFLATION SYMBOL

Inflation symbol disappears and inflation becomes slow.

PULSE SYMBOL is displayed as pulse is detected.


PULSE SYMBOL

Air is released from the cuff and measured values are displayed as measurement completes.

Refer to page 3 for [].

Refer to page 4 for [].

Refer to page 12 for "Err" display.

4 SELECT MEMORY BANK TO SAVE THE MEASUREMENT RESULT.

MEMORY SYMBOL is displayed below MEMORY KEY.

Touching MEMORY KEY will switch the bank.

Readings are automatically saved in the bank, either MEMORY [1] or [2], selected when the monitor is turned off after measurement. Make sure to confirm the bank selected after a measurement.

Refer to page 4 for memory function.

5 TURN OFF THE MONITOR.

Touch "START/STOP" KEY.

The monitor will be turned off. Even if you do not turn off the monitor, it turns off automatically after 3 minutes.

Do not execute repeated measurements for congestion of blood could result in false measurement. Let your arm rest for at least 5 minutes.

VIEWING THE SAVED READINGS (MEMORY FUNCTION)

Turn off the monitor once with "START/STOP" KEY when the results are displayed after a measurement.

Refer to page 4 for memory function.


To delete readings, refer to page 11.

1 TOUCH MEMORY KEY.

Touch MEMORY KEY [1] to review readings saved in the memory bank [1] and MEMORY KEY [2] to review those in the bank [2].

The average of the saved readings are displayed.

Note: Average will not be displayed unless there are two or more readings saved.

A : Average
NUMBER OF SAVED RESULTS


2 TOUCH THE SAME MEMORY KEY AGAIN.

The latest reading saved in the selected bank is displayed.

The display at the top shows memory number, measurement date and time alternately. The latest reading is displayed as memory "No. 1".

The results are recalled from the latest to the oldest with operation of MEMORY KEY.

The bigger memory number indicates the older number.

SWITCHING MEMORY BANKS

Touching MEMORY KEY [2] at display of reading in the bank [1] will switch the display to reading in bank [2] and touching MEMORY KEY [1] again will return to the display of result in the bank [1].

3 TURN OFF THE MONITOR.

TROUBLESHOOTING

??? Blood pressure is extremely high or low.

Measurement is not taken with correct posture or the cuff is wrapped incorrectly.
Reconfirm the measurement procedures.

The reading was effected from moving or talking. **Remain still and quiet during measurement.**

Measurement was taken right after exercise or other activities influencing the reading. **Measure again after resting for more than 5 minutes. See TIPS ON TAKING YOUR BLOOD PRESSURE.**

Touch "START/STOP" KEY.

Even if you do not turn off the monitor, it displays the result for approximately 30 seconds and will turn off automatically.

DELETING SAVED READINGS

Deleting individual readings

Touch the memory key to display the latest set of readings. Then scroll through the readings by repeatedly touching the memory key until you reach the reading to be deleted. Then hold down memory key for 4 seconds until display flashes and continue to hold down key for about 4 seconds until reading is deleted and the display is cleared [- - -].

Deleting all the readings in a memory bank.

Select and display the average readings for the memory bank to be deleted. Then delete by holding down the memory key continuously until the display is cleared [- - -] as described in the sentence above for the deletion of individual readings.

??? Measured values vary all the time.

Mental and physical condition or measurement condition is influencing blood pressure. **Make measurement under the same condition.**

??? Measured value is different from that taken at a hospital.

Mental condition such as nervousness at a hospital or relaxed feeling at home is influencing blood pressure. **Take records of blood pressure taken at home and consult with your doctor.**

Blood pressure is changing 24 hours a day in relation with mental and physical conditions. Eating, drinking, smoking, exercising or bathing will alter blood pressure. Blood pressure will be also influenced by mental condition such as feeling tensed or relieved and health condition. It is important to know the tendency of blood pressure by making measurement regularly. Make your own rules including conditions and time to make measurement and take blood pressure everyday.

[▲] symbol is displayed again in the course of inflation and deflation and inflation of the cuff are repeated.

When the initial pressure was too high for blood pressure to be taken or if you move or strain your arm or hand, the cuff is inflated again. Deflation and inflation of the cuff are repeated until blood pressure is taken. This does not indicate any problem or malfunction of the monitor. **Do not move or speak during measurement.**

[Err 300] is displayed.

Over pressure: blood pressure could not be taken due to moving or speaking although the cuff was inflated to the maximum pressure. **Do not move or speak during measurement.**

[Err -1] is displayed.

Blood pressure could not be taken due to moving or speaking. **Do not move or speak during measurement.**

[Err -2] is displayed.

The cuff is not securely connected to the monitor. Secure the air connector.

The cuff is not wrapped appropriately. **Reconfirm how to put the cuff on your arm.**

[?] is displayed.

Batteries are weak. Replace all batteries with new ones.

Nothing is displayed.

Batteries are depleted. **Replace all batteries with new ones.**

Batteries are inserted incorrectly. **Reinsert the batteries correctly.**

Battery terminals are not clean. **Clean the terminals with dry cloth.**

AC adaptor is not connected. **Connect AC adaptor.**

Touching the keys with tip of finger or nail is touching the key. **Touch the key with the plane of your finger.**

Clock is not displayed.

The clock is inactivated. NOTE: Clock is inactivated when batteries and/or AC adaptor are/is removed. **Adjust date and time and activate the clock.**

Measurement date and time are displayed with [---] and [-:-].


The clock is not activated. **Adjust date and time and activate the clock.**

Measurement was taken before the clock was activated. Measurement date and time cannot be saved without clock being activated.


[00] is displayed.

"START/STOP" KEY was being touched when batteries are being installed or when AC adaptor is being plugged in. **Turn off the monitor once with "START/STOP" KEY.**

Cuff end is loosened.

Put the cuff with its hook and loop fastener side down. Put the cuff end through the metal ring. Fold back the cuff at the metal ring.

If you cannot get correct measurement with the methods above, contact your dealer. Do

TECHNICAL DESCRIPTION

DSK-1011 complies with the EMC, electromagnetic compatibility, standard, IEC60601-1-2. Refer to the following tables for specific information regarding compliance to the standard.

DSK-1011, as a medical electrical equipment, needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

- Portable and mobile RF communications equipments can affect the device.
- The use of accessories other than those specified in this manual may result in increased emissions or decreased immunity of the device.
- DSK-1011 should not be used adjacent to or stacked with other equipment.

Table 201 - Guidance and manufacturer's declaration - electromagnetic emissions

DSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1011 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	DSK-1011 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	DSK-1011 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 202 - Guidance and manufacturer's declaration - electromagnetic immunity

DSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1011 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _r (>95% dip in U _r) for 0.5 cycle 40% U _r (60% dip in U _r) for 5 cycles 70% U _r (30% dip in U _r) for 25 cycles <5% U _r (>95% dip in U _r) for 5 sec	<5% U _r (>95% dip in U _r) for 0.5 cycle 40% U _r (60% dip in U _r) for 5 cycles 70% U _r (30% dip in U _r) for 25 cycles <5% U _r (>95% dip in U _r) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of DSK-1011 requires continued operation during power mains interruptions, it is recommended that DSK-1011 is be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE	U _r is the a.c. mains voltage prior to application of the test level.		

not disassemble or tamper with the internal mechanism.

Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity

DSK-1011 is intended for use in the electromagnetic environment specified below. The customer or the user of DSK-1011 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of DSK-1011, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=[3.5/V]^P$
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m	$d=[3.5/E]^P$, 80MHz-800MHz $d=[7/E]^P$, 800MHz-2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ^a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which DSK-1011 is used exceeds the applicable RF compliance level above, DSK-1011 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating DSK-1011.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and DSK-1011

DSK-1011 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of DSK-1011 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and DSK-1011 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz to 80 MHz, $d=[3.5/V]^P$	80 MHz to 800 MHz, $d=[3.5/E]^P$	800 MHz to 2.5 GHz, $d=[7/V]^P$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

This equipment is guaranteed for the period of 2 years after the date of purchase against manufacturing defects when returned along with the proof of date of purchase to the dealer from whom the purchase was made. During this period, the unit will be repaired or replaced free of charge if the failure is attributable to faulty design or manufacture. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual. Please contact your dealer for additional warranty provisions which may remain effective after the manufacturer's warranty period has expired.

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EC-Representative:
Nissei Healthcare (UK) Ltd.
Henfield, BN5 9SJ UK
web site <http://www.nisseihealthcare.com>

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NIHON SEIMITSU SOKKI CO., LTD.