Spengler

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USER'S MANUAL

MANUAL

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USER'S MANUAL

EN

Thank you for purchasing a MAXI+3 large dial sphygmomanometer. When you choose a SPENGLER instrument, you are choosing experience and excellence that has been recognized since 1907.

INSTRUCTIONS FOR USE

Please read the leaflet before use and keep it for future reference.

INTENDED MEDICAL USE

The MAXI +3 is a medical diagnosis professional device (used by health professionals trained in its use) that measures the pressure exercised by blood on arteries by the auscultatory method (a stethoscope is therefore required). It accurately measures maximum systolic blood pressure or SBP (resulting from the left ventricle contracting) and minimum diastolic blood pressure or DBP (resulting from ventricles relaxing).

WHAT DOES "BLOOD PRESSURE" MEAN?

Blood pressure is the pressure exercised by blood flow on the walls of arteries.The highest pressure, produced in the arteries by each heart beat, is called "systolic pressure" whereas the lowest is "diastolic pressure".

Blood pressure fluctuates several times a day and can also vary with stress, nutritional factors, physical activity, medicines or diseases.

WHAT IS "NORMAL" BLOOD PRESSURE?

According to the WHO (World Health Organization), normal blood pressure is systolic pressure less than 130mmHg (millimetres of mercury) and diastolic pressure below 90mmHg. However, blood pressure can vary depending on the individual. The WHO classifications table can be found below.

Blood pressure classification in adults (mmHg)

	Systolic blood pressure (mmHg)	And/or	Diastolic blood pressure (mmHg)
Optimum	< 120	and	< 80
Normal	120-129	and/or	80-84
Upper normal	130-139	and/or	85-89
Grade 1 HBP	140-159	and/or	90-99
Grade 2 HBP	160-179	and/or	100-109
Grade 3 HBP	≥ 180	and/or	≥ 110

ESH 2007 Recommendations.

WARNINGS

The warnings in this document highlight conditions or practices which could cause injury, disease or potentially death of a patient.

⚠ **Warning:** Possible measurement error.Only use the tensiometer armbands and MAXI+3 accessories from SPENGLER to avoid measurement errors.

⚠ **Warning:** Risk of inaccurate measurement. Before use, check that all the connection points, the cuff and tube are completely airtight. Excessive leakage can affect measurement.

⚠ **Warning:** No changes can be made to this equipment as this could affect measurements.

CHOOSE AN CUFF THAT SUITS THE PATIENT'S MORPHOLOGY

SPENGLER is taking on the fight by informing health professionals in using an cuff that suits a patient's morphology, as the size of the cuff impacts the quality of measurements:

- A too narrow cuff over-estimates blood pressure.
- A too large cuff could under-estimate blood pressure.

Note:

The index marker showing the arterial index on the cuff must be within the range shown on Maxi+3 cuffs.

When the cuff is positioned on the patient's arm, the "index" marker shown on the cuff must be within the range shown on Maxi+3 cuffs. If the "index" marker or the edge of the marker does not reach this range, we recommend that you see the table below in order to adapt the cuff to a patient's morphology.

Available sizes of armbands for MAXI+3 SPHYGMOMANOMETER

Size of cuff	Arm circumference
S	18,4 cm – 29,7 cm
M	25,4 cm – 40,6cm
L	33,3 cm - 51 cm

USING THE MAXI+3 LARGE DIAL ANEROID SPHYGMOMANOMETER

Steps

- Leave the patient to rest for 5 minutes
- Place the patient in a comfortable position, with legs uncrossed, and the back and arm elevated, ideally without moving or speaking
- Check that the marker index for arterial index is within the range shown on the Maxi+3 cuffs
- Keep the cuff at heart height during measurement
- Keep the pneumatic circuit closed (check the decompression valve)
- Inflate the cuff for an initial SBP estimation (when the blood circulates again, it causes a noticeable beating sound on the stethoscope) until radial pulse disappears (so-called Korotkoff sounds)
- Inflate again to 30mmHg above the estimated SBP
- Deflate by opening the valve at a speed of 2mmHg per second and simultaneously read the blood pressure on the manometer
- Read the SBP for phase I of Korotkoff sounds
- Read the DBP for phase V of Korotkoff sounds (complete disappearance of sounds noted at the stethoscope)
- Measure the heart rate
- Continue decompression until complete purge of the cuff

⚠ **Warning:** A blood flow stoppage carries risks.

Recommandation:

- Use K5* for adult auscultation
- Use K5* for auscultation of pregnant patients, so long as the sounds are noticeable when the cuff is deflated, in which case the K4* can be used
- Use K4* for auscultation of children between 3 and 12 years.

** See the table describing the phases of Korotkoff sounds.*

KOROTKOFF SOUNDS

Environmental specifications:

Operating temperature: +10°C to +40°C

Operating relative humidity: 20% to 85%

Storage temperature: -20°C to +70°C

Relative humidity at storage: 20% to 85%

MAINTENANCE AND MAINTENANCE VERIFICATION

Comply with current regional regulations for the disposal of manual tensiometers or accessories.

Metrological inspection by an approved body

Check accuracy and watertightness:

- At least once every 2 years (manometer)
- At least once per year (cuffs, tube, bulb)
- After each maintenance operation or repair or replacement of removable parts.

Note: In accordance with international and national metrological standards, the operator or the body responsible for maintenance must ensure traceability of benchmark manometer used for calibration.

Removable parts & accessories:

Please contact your approved retailer to purchase the parts and accessories below:

- Replacement switch,
- Tube,
- Armbands (all sizes).

CLEANING AND DISINFECTION:

⚠ **Warning:** Never immerse the manometer in liquid or in a decontamination device. Use a neutral detergent to clean the surface of product. Should disinfection be necessary, use a cloth soaked in 75% alcohol.

⚠ **Warning:** Avoid placing your sphygmomanometer too close to heat, excessive cold, solvents or greases. Avoid long-term contact with sweat, as the PVC could easily stiffen (switch and tube). Keep hands dry.

TECHNICAL SPECIFICATIONS:

Maxi+3 Manometer

Accuracy: +/-3mmHg

Operating range:0 to 300mmHg

Displayed range:0 to 300mmHg

Indication type: radial scale

Scale grading: increments of 2 mmHg

Pressure generation: switch

Pressure reduction: air release valve

Case:ABS

Mechanism: copper-beryllium alloy

Switch:PVC

Tube: PVC.

Cuff: Nylon.

PRECAUTIONS FOR USE

Avoid impacts, particularly on the dial glass.

Do not expose to heat, sun or humidity, dust.

WARRANTY

The MAXI+3 sphygmomanometer is guaranteed against any part and labour or operating fault, in accordance with the manufacturer's specifications, under normal maintenance and usage conditions.The 2-year guarantee for the tensiometer (excluding armband and tube) starts from the date of purchase from an approved retailer. SPENGLER's obligations are limited during the guarantee period to the repair or replacement of components deemed to be defective by SPENGLER.

Guarantee for Maxi+3 accessories:

Tube:1 year guarantee

Inflation switch:1 year guarantee

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STANDARD

Maxi+3 is compliant with the European reference standard NF EN ISO 81060-1.

QUALITY ASSURANCE

NF EN ISO 13485 :2016

CE MARKING

CE certification aneroid sphygmomanometer provided by TUV SUD, notified body no. 0123. Class Im device (measurement of blood pressure by auscultatory method) in accordance with European directive 93/42/CEE regarding medical devices.

Symbol	Description
	Manufacturer
	See the directions for use
	Product reference
	Batch number
	Protect from humidity
	Protect from light and sun
	Storage temperature range
	Dispose of used product at a collection point for recycling in accordance with local regulations
	This product complies with directive 93/42/CEE regarding medical devices.

DESIGNATION	REFERENCE
Maxi+3 Wall mounted with adult cuff (M)	522 700
Maxi+3 with stand with adult cuff (M)	522 705
Maxi+3 Wall mounted with multicuff (L, M, S)	522 707
Maxi+3 with stand with multicuff (L, M, S)	522 708

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