

# Solutions for Fragile Skin



3M™ Red Dot™ Soft Cloth Monitoring Electrode 2238



3M™ Red Dot™ Monitoring Electrode with 3M™ Micropore™ Tape Backing 2239



3M™ Red Dot™ ECG Monitoring Electrode 2268-3



3M™ Red Dot™ Monitoring Electrode with 3M™ Micropore™ Tape Backing 2249-50 front and back with blue abradant



3M™ Red Dot™ Soft Cloth Monitoring Electrode 2255-50 front and back with blue abradant

In today's patient population, fragile skin is often a concern whether due to age, medication or repeated removal of adhesive products. It's a delicate balance to protect the patient's skin and maintain good adhesion to obtain acceptable EKG trace quality.

**Conformable backings to follow skin contours.** 3M™ Micropore™ and soft cloth tape backings are low profile to minimize catching of medical devices and fabrics. This reduces the risk of early lift which often necessitates electrode removal and replacement.

**Gentle adhesives.** Pressure Sensitive Adhesives (PSAs) need both pressure and warmth (from the patient's body) to achieve optimum adhesion. 3M has formulated these electrodes PSAs to be gentler to the skin.

**Small size option.** The 2268-3 electrode is a compact addition to our fragile skin portfolio, well-suited for both pediatric and adult patients. Gentle, conductive adhesive over the entire electrode may improve trace quality. Low profile increases conformability and comfort. Lift tab makes electrode removal easy.

## Suggested Applications

	2238	2239	2249-50	2255-50	2268-3
Cath Lab	●	●	●	●	●
Echo-Cardio	●	●	●	●	●
ICU	●	●	●	●	●



To learn more about 3M™ Patient Monitoring Solutions,  
call the **3M Health Care Helpline** at **1-800-228-3957**  
or visit **[www.3M.com/patientmonitoring](http://www.3M.com/patientmonitoring)**.



**3M Health Care**  
**Infection Prevention Division**  
3M Center, Building 275-4E-01  
St. Paul, MN 55144-1000  
USA  
1-800-228-3957  
[www.3M.com/patientmonitoring](http://www.3M.com/patientmonitoring)

Please recycle. Printed in U.S.A.  
© 3M 2014. All rights reserved.  
70-2008-9621-8 HB9710

Philips and IntelliVue are  
trademarks of Koninklijke  
Philips Electronics N.V.  
Limited Liability Company.

GE is a trademark of General  
Electric Company.

GE Healthcare, a division of  
General Electric Company.

Spacelabs Healthcare is a  
wholly owned subsidiary of  
OSI Systems, Inc.

Nihon Kohden is a trademark  
of Nihon Kohden Corporation.

Datascope is a registered  
trademark of Datascope  
Corporation.

Siemens/Dräger is  
trademark of Siemens  
Aktiengesellschaft.

3M, Red Dot and Micropore  
are trademarks of 3M.