

CONFIGURATION PC REQUISE POUR LOGICIELS CARDIOLINE

CONFIGURATIONS HARDWARE REQUISES

TABLETTE PC		
Logiciels Cardioline	Specifications Techniques	Minimum requis
TouchECG Glasgow interpretation WebUploader DICOM	• Processeur • Affichage • RAM/HDD • OS • Input / Output • Communication	Processeur Intel Atom Quad core 1.6 GHz ou plus Touch-screen, 800x600 minimum 4 GB / 64 GB Windows 7, 8.1 or 10 with netframework 3.5 (touch) USB Bluetooth (Win BT stack) /WiFi e/o 3G

All-in-One PC		
Logiciels Cardioline	Caractéristiques	Minimum requis
TouchECG Glasgow interpretation Cubesuite CubeStress WebUploader / WebApp pers ou Private Pneumos500 PC App DICOM	• Processeur • Affichage • RAM/HDD • OS • Accessoires • Input / Output • Communication	Processeur Intel i3 ou au-delà 1366 x 768 (écran tactile souhaitable) 4 GB / 64 GB Windows 7, 8.1 or 10 with netframework 3.5 (touch) Clavier et souris, DVD, micro et HP USB / RS232 (application stress test) Bluetooth (Win BT stack) /WiFi/LAN